

# CRB High Throughput Lab in Riverside For HLB Diagnostics




Cynthia LeVesque  
Laboratory Director  
Citrus Research Board



# Lab Location

- Suite 500, 1201 Research Park Drive  
– corner of Columbia and Research Park Drive
- UC Riverside campus 5 minutes away



UC  
Riverside







Laboratory space was a single large room.

Has been divided into:


- Sample processing
- DNA extraction / PCR set up rooms






Small PCR instrument room.

Has been:

- Isolated by a hallway from other lab areas
- Entered only by lab personnel



UR = unisex restroom  
H = hood  
S = sink  
A = autoclave

FH = flow hood  
 = new benches  
 = new walls  
 = new doors

# Lab renovation: 1<sup>st</sup> steps



4 Preparing for "mud" application



3 Dry wall installation



1 Framing the hallway



2 Framing the sample lab



# Lab renovation: relocation of the fume hood and installation of the ice maker



3 Ice maker installed and fittings for RO water system



1 Hood moved for larger sink and ice maker



2 Installation of new counter next to sink

# Lab renovation: Installation of counters and cabinets





# Sample Flow

---

## Reception

- Download from data logger / data entry from data forms

Associates samples with:

Tree bar code

GPS coordinates

Field photographs

- Scan data forms

Electronically captures paper forms

- Photograph samples

Accumulation of symptom library associated with test results



- Bar code sample

Associates sample with tree

Associates RT-PCR results with sample

# Sample Flow

## Sample Prep Room

- Dissect plant tissue into 1 ml 96 well plate  
200 mg of leaf midvein and petiole
- Autoclave remaining plant tissue  
Eliminates possibility of contamination from plant material
- Freeze dry plant samples  
Stabilizes sample at room temperature
- Collect live insect samples  
Open sample with insects inside insect cages  
Remove psyllid sample into EtOH



- Remove insects from sticky traps using HistoClear
- Freeze all remaining tissue to kill insects



# Sample Flow

---

## DNA / PCR Prep Room

- Pulverize freeze dried plant tissue  
96 well format bead dispenser and beater
- Extract DNA  
Initially will use USDA certified procedure using DNeasy Plant Kit (Quiagen, Valencia, CA)  
Evaluating 96 well format robot for DNA extraction, normalization, and RT-PCR set up
- Set up RT-PCR reactions in high performance laminar flow hood

## Instrument Room

- Run RT-PCR on ABI 7500 Fast system using USDA certified primers and parameters  
Run samples will not re-enter the DNA / PCR Prep Room  
Samples will be bagged for biohazard waste disposal

# Currently

---

- Construction completed, city certifications in progress
- All equipment in place by the end of June
- ABI 7500 Fast installation scheduled for June 26<sup>th</sup>
- Ready to perform DNA extraction of citrus and insect samples
- Running RT-PCR reactions by mid July
- Will seek USDA certification by the end of the summer